

Product Datasheet

Mouse Glucose-6-phosphatase ELISA Kit (orb1212216)

Description	Mouse Glucose-6-phosphatase, G6PC ELISA Kit
Reactivity	Mouse
Range	0.5-150U/mL
Tested Applications	ELISA
Target	G6PC
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	40,473
Uniprot ID	P35576
Solubility (25°C)	160U/mL
Sensitivity	0.242U/mL
Entrez	14377
Expiration Date	Please enquire.