

Product Datasheet

Mouse Glucose Regulated Protein 78 ELISA Kit (orb1212217)

Description	Mouse Glucose Regulated Protein 78, GRP78 ELISA Kit
Reactivity	Mouse
Range	0.2-60ng/mL
Tested Applications	ELISA
Target	HSPA5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	72,422
Uniprot ID	P20029
Solubility (25°C)	64ng/mL
Sensitivity	0.121ng/mL
Entrez	14828
Expiration Date	Please enquire.