

Product Datasheet

Mouse Gdnf Family Receptor Alpha-like ELISA Kit (orb1212219)

Description	Mouse Gdnf Family Receptor Alpha-like, GFRAL ELISA Kit
Reactivity	Mouse
Range	15-3000ng/L
Tested Applications	ELISA
Target	GFRAL
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	43,964
Uniprot ID	Q6SJE0
Solubility (25°C)	3200ng/L
Sensitivity	7.85ng/L
Entrez	404194
Expiration Date	Please enquire.