

Product Datasheet

Mouse Galanin-like Peptide ELISA Kit (orb1212226)

Description	Mouse Galanin-like Peptide, GALP ELISA Kit
Reactivity	Mouse
Range	15-3000ng/L
Tested Applications	ELISA
Target	GALP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	12,773
Uniprot ID	Q810H5
Solubility (25°C)	3200ng/L
Sensitivity	7.47ng/L
Entrez	232836
Expiration Date	Please enquire.