

Product Datasheet

Mouse Fibulin-1 ELISA Kit (orb1212231)

Description	Mouse Fibulin-1, FBLN1 ELISA Kit
Reactivity	Mouse
Range	1-400ng/mL
Tested Applications	ELISA
Target	FBLN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	78,033
Uniprot ID	Q08879
Solubility (25°C)	480ng/mL
Sensitivity	0.56ng/mL
Entrez	14114
Expiration Date	Please enquire.