

Product Datasheet

Mouse Fibroblast Growth Factor 8 ELISA Kit (orb1212232)

Description	Mouse Fibroblast Growth Factor 8, FGF8 ELISA Kit
Reactivity	Mouse
Range	3-900ng/L
Tested Applications	ELISA
Target	FGF8
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	30,420
Uniprot ID	P37237
Solubility (25°C)	960ng/L
Sensitivity	1.46ng/L
Entrez	14179
Expiration Date	Please enquire.