

Product Datasheet

Mouse Fibroblast Growth Factor 5 ELISA Kit (orb1212233)

Description	Mouse Fibroblast Growth Factor 5, FGF5 ELISA Kit
Reactivity	Mouse
Range	20-4500ng/L
Tested Applications	ELISA
Target	FGF5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	29,103
Uniprot ID	P15656
Solubility (25°C)	4800ng/L
Sensitivity	9.29ng/L
Entrez	14176
Expiration Date	Please enquire.