

Product Datasheet

Mouse Epidermal Growth Factor-like Protein 6 ELISA Kit (orb1212238)

Description	Mouse Epidermal Growth Factor-like Protein 6, EGFL6 ELISA Kit
Reactivity	Mouse
Range	50-10000ng/L
Tested Applications	ELISA
Target	EGFL6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	61,520
Uniprot ID	Q9JJZ5
Solubility (25°C)	12800ng/L
Sensitivity	24.63ng/L
Entrez	54156
Expiration Date	Please enquire.