

Product Datasheet

Mouse Elongation Factor 2 ELISA Kit (orb1212242)

Description	Mouse Elongation Factor 2, EEF2 ELISA Kit
Reactivity	Mouse
Range	20-6000ng/L
Tested Applications	ELISA
Target	EEF2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	95,314
Uniprot ID	P58252
Solubility (25°C)	6400ng/L
Sensitivity	9.84ng/L
Entrez	13629
Expiration Date	Please enquire.