

## Product Datasheet

# Mouse EGF-containing Fibulin-like Extracellular Matrix Protein 1 ELISA Kit (orb1212246)

|                            |                                                                                    |
|----------------------------|------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse EGF-containing Fibulin-like Extracellular Matrix Protein 1, EFEMP1 ELISA Kit |
| <b>Reactivity</b>          | Mouse                                                                              |
| <b>Range</b>               | 20-4500ng/L                                                                        |
| <b>Tested Applications</b> | ELISA                                                                              |
| <b>Target</b>              | EFEMP1                                                                             |
| <b>Storage</b>             | 2-8°C                                                                              |
| <b>Note</b>                | For research use only                                                              |
| <b>Application notes</b>   | Sandwich                                                                           |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates                                             |
| <b>MW</b>                  | 54,953                                                                             |
| <b>Uniprot ID</b>          | <b>Q8BPB5</b>                                                                      |
| <b>Solubility (25°C)</b>   | 4800ng/L                                                                           |
| <b>Sensitivity</b>         | 9.95ng/L                                                                           |
| <b>Entrez</b>              | <b>216616</b>                                                                      |
| <b>Expiration Date</b>     | Please enquire.                                                                    |