

Product Datasheet

Mouse Drebrin ELISA Kit (orb1212252)

Description	Mouse Drebrin, DBN1 ELISA Kit
Reactivity	Mouse
Range	2-600ng/L
Tested Applications	ELISA
Target	DBN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	77,287
Uniprot ID	Q9QXS6
Solubility (25°C)	640ng/L
Sensitivity	0.97ng/L
Entrez	56320
Expiration Date	Please enquire.