

Product Datasheet

Mouse Disintegrin and Metalloproteinase Domain-containing Protein 12 ELISA Kit (orb1212255)

Description	Mouse Disintegrin and Metalloproteinase Domain-containing Protein 12, ADAM12 ELISA Kit
Reactivity	Mouse
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	ADAM12
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	98,504
Uniprot ID	Q61824
Solubility (25°C)	24ng/mL
Sensitivity	0.024ng/mL
Entrez	11489
Expiration Date	Please enquire.