

Product Datasheet

Mouse Polycomb Complex Protein Bmi-1 ELISA Kit (orb1212258)

Description	Mouse Polycomb Complex Protein Bmi-1, BMI-1 ELISA Kit
Reactivity	Mouse
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	BMI1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	36,708
Uniprot ID	P25916
Solubility (25°C)	32ng/ml
Sensitivity	0.022ng/ml
Entrez	12151
Expiration Date	Please enquire.