

Product Datasheet

Mouse Apolipoprotein A1 ELISA Kit (orb1212262)

Description	Mouse Apolipoprotein A1, APO-A1 ELISA Kit
Reactivity	Mouse
Range	20-4500ng/mL
Tested Applications	ELISA
Target	APOA1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	30,616
Uniprot ID	Q00623
Solubility (25°C)	4800ng/mL
Sensitivity	9.01ng/mL
Entrez	11806
Expiration Date	Please enquire.