

Product Datasheet

Mouse Chordin-like Protein 1 ELISA Kit (orb1212290)

Description	Mouse Chordin-like Protein 1, CHRD1 ELISA Kit
Reactivity	Mouse
Range	10-2000ng/L
Tested Applications	ELISA
Target	CHRD1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	50,732
Uniprot ID	Q920C1
Solubility (25°C)	2400ng/L
Sensitivity	5.15ng/L
Entrez	83453
Expiration Date	Please enquire.