

Product Datasheet

Mouse C-C Chemokine Receptor Type 2 ELISA Kit (orb1212298)

Description	Mouse C-C Chemokine Receptor Type 2, CCR2 ELISA Kit
Reactivity	Mouse
Range	10-2000ng/L
Tested Applications	ELISA
Target	CCR2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	42,783
Uniprot ID	P51683
Solubility (25°C)	2400ng/L
Sensitivity	5.49ng/L
Entrez	12772
Expiration Date	Please enquire.