

Product Datasheet

Mouse Caveolin-1 ELISA Kit (orb1212300)

Description	Mouse Caveolin-1, CAV1 ELISA Kit
Reactivity	Mouse
Range	50-10000ng/L
Tested Applications	ELISA
Target	CAV1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,539
Uniprot ID	P49817
Solubility (25°C)	12800ng/L
Sensitivity	23.57ng/L
Entrez	12389
Expiration Date	Please enquire.