

Product Datasheet

Mouse Bone Marrow Proteoglycan ELISA Kit (orb1212307)

Description	Mouse Bone Marrow Proteoglycan, PRG2 ELISA Kit
Reactivity	Mouse
Range	1-400ng/mL
Tested Applications	ELISA
Target	PRG2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	24,255
Uniprot ID	Q61878
Solubility (25°C)	480ng/mL
Sensitivity	0.53ng/mL
Entrez	19074
Expiration Date	Please enquire.