

Product Datasheet

Mouse Beta-2 Adrenergic Receptor ELISA Kit (orb1212310)

Description	Mouse Beta-2 Adrenergic Receptor, ADRB2 ELISA Kit
Reactivity	Mouse
Range	10-2000ng/L
Tested Applications	ELISA
Target	ADRB2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	46,998
Uniprot ID	P18762
Solubility (25°C)	2400ng/L
Sensitivity	5.29ng/L
Entrez	11555
Expiration Date	Please enquire.