

Product Datasheet

Mouse Beta-1 Adrenergic Receptor ELISA Kit (orb1212311)

Description	Mouse Beta-1 Adrenergic Receptor, ADRB1 ELISA Kit
Reactivity	Mouse
Range	20-6000ng/L
Tested Applications	ELISA
Target	ADRB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	50,494
Uniprot ID	P34971
Solubility (25°C)	6400ng/L
Sensitivity	9.75ng/L
Entrez	11554
Expiration Date	Please enquire.