

Product Datasheet

Mouse Nt-3 Growth Factors Receptor ELISA Kit (orb1212312)

Description	Mouse Nt-3 Growth Factors Receptor, NTRK2 ELISA Kit
Reactivity	Mouse
Range	20-4500ng/ml
Tested Applications	ELISA
Target	NTRK2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	92,133
Uniprot ID	P15209
Solubility (25°C)	4800ng/ml
Sensitivity	9.23ng/ml
Entrez	18212
Expiration Date	Please enquire.