

Product Datasheet

Mouse Apolipoprotein A-4 ELISA Kit (orb1212317)

Description	Mouse Apolipoprotein A-4, APOA4 ELISA Kit
Reactivity	Mouse
Range	7-1500ng/mL
Tested Applications	ELISA
Target	APOA4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	45,029
Uniprot ID	P06728
Solubility (25°C)	1600ng/mL
Sensitivity	3.13ng/mL
Entrez	11808
Expiration Date	Please enquire.