

Product Datasheet

Mouse Allograft Inflammatory Factor 1 ELISA Kit (orb1212328)

Description	Mouse Allograft Inflammatory Factor 1, AIF1 ELISA Kit
Reactivity	Mouse
Range	15-3000ng/L
Tested Applications	ELISA
Target	AIF1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	16,911
Uniprot ID	070200
Solubility (25°C)	3200ng/L
Sensitivity	7.75ng/L
Entrez	11629
Expiration Date	Please enquire.