

Product Datasheet

Mouse Zinc Finger Protein Gli1 ELISA Kit (orb1212334)

Description	Mouse Zinc Finger Protein Gli1, GLI1 ELISA Kit
Reactivity	Mouse
Range	0.03-9ng/ml
Tested Applications	ELISA
Target	GLI1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	118,560
Uniprot ID	P47806
Solubility (25°C)	9.6ng/ml
Sensitivity	0.016ng/ml
Entrez	14632
Expiration Date	Please enquire.