

Product Datasheet

Mouse Cyclin Dependent Kinase Inhibitor 3 ELISA Kit (orb1212345)

Description	Mouse Cyclin Dependent Kinase Inhibitor 3, CDKN3 ELISA Kit
Reactivity	Mouse
Range	15-3000ng/L
Tested Applications	ELISA
Target	CDKN3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	23,793
Uniprot ID	Q810P3
Solubility (25°C)	3200ng/L
Sensitivity	7.86ng/L
Entrez	72391
Expiration Date	Please enquire.