

Product Datasheet

Mouse Fibroblast Growth Factor 9 ELISA Kit (orb1212354)

Description	Mouse Fibroblast Growth Factor 9, FGF9 ELISA Kit
Reactivity	Mouse
Range	1-400ng/L
Tested Applications	ELISA
Target	FGF9
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	23,414
Uniprot ID	P54130
Solubility (25°C)	480ng/L
Sensitivity	0.48ng/L
Entrez	14180
Expiration Date	Please enquire.