

Product Datasheet

Mouse Aryl Hydrocarbon Receptor ELISA Kit (orb1212356)

Description	Mouse Aryl Hydrocarbon Receptor, AHR ELISA Kit
Reactivity	Mouse
Range	6-1500ng/L
Tested Applications	ELISA
Target	AHR
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	95,017
Uniprot ID	P30561
Solubility (25°C)	1600ng/L
Sensitivity	3.22ng/L
Entrez	11622
Expiration Date	Please enquire.