

Product Datasheet

Mouse Cluster of Differentiation 4 ELISA Kit (orb1212359)

Description	Mouse Cluster of Differentiation 4, CD4 ELISA Kit
Reactivity	Mouse
Range	3-900ng/ml
Tested Applications	ELISA
Target	NEUROD4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	37,133
Uniprot ID	009105
Solubility (25°C)	960ng/ml
Sensitivity	0.16ng/ml
Entrez	11923
Expiration Date	Please enquire.