

Product Datasheet

Mouse Fibroblast Growth Factor 2 ELISA Kit (orb1212360)

Description	Mouse Fibroblast Growth Factor 2, FGF-2 ELISA Kit
Reactivity	Mouse
Range	3-900ng/L
Tested Applications	ELISA
Target	FGF2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	17,153
Uniprot ID	P15655
Solubility (25°C)	960ng/L
Sensitivity	0.14ng/L
Entrez	14173
Expiration Date	Please enquire.