

Product Datasheet

Mouse Solute Carrier Family 12 Member 2 ELISA Kit (orb1212361)

Description	Mouse Solute Carrier Family 12 Member 2, SLC12A2 ELISA Kit
Reactivity	Mouse
Range	5-1800ng/L
Tested Applications	ELISA
Target	SLC12A2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	131,034
Uniprot ID	P55012
Solubility (25°C)	1920ng/L
Sensitivity	2.38ng/L
Entrez	20496
Expiration Date	Please enquire.