

Product Datasheet

Mouse Microtubule-associated Protein 1 Light Chain 3 Alpha ELISA Kit (orb1212362)

Description	Mouse Microtubule-associated Protein 1 Light Chain 3 Alpha, MAP1LC3A ELISA Kit
Reactivity	Mouse
Range	10-4000ng/L
Tested Applications	ELISA
Target	MAP1LC3A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	14,272
Uniprot ID	Q91VR7
Solubility (25°C)	4800ng/L
Sensitivity	5.24ng/L
Entrez	66734
Expiration Date	Please enquire.