

Product Datasheet

Mouse Dieptidyl Peptidase 4 ELISA Kit (orb1212363)

Description	Mouse Dieptidyl Peptidase 4, DPP-4 ELISA Kit
Reactivity	Mouse
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	DPP4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	87,437
Uniprot ID	P28843
Solubility (25°C)	40ng/ml
Sensitivity	0.03ng/ml
Entrez	13482
Expiration Date	Please enquire.