

Product Datasheet

Mouse Leukotriene A4 Hydrolase ELISA Kit (orb1212364)

Description	Mouse Leukotriene A4 Hydrolase, LTA4H ELISA Kit
Reactivity	Mouse
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	LTA4H
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	69,051
Uniprot ID	P24527
Solubility (25°C)	16ng/ml
Sensitivity	0.027ng/ml
Entrez	16993
Expiration Date	Please enquire.