

Product Datasheet

Mouse Sestrin 2 ELISA Kit (orb1212373)

Description	Mouse Sestrin 2, SESN2 ELISA Kit
Reactivity	Mouse
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	SESN2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	54,377
Uniprot ID	P58043
Solubility (25°C)	32ng/ml
Sensitivity	0.027ng/ml
Entrez	230784
Expiration Date	Please enquire.