

Product Datasheet

Mouse Tyrosinase ELISA Kit (orb1212375)

Description	Mouse Tyrosinase, TYR ELISA Kit
Reactivity	Mouse
Range	20-4500ng/L
Tested Applications	ELISA
Target	TYR
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	60,606
Uniprot ID	P11344
Solubility (25°C)	4800ng/L
Sensitivity	9.24ng/L
Entrez	22173
Expiration Date	Please enquire.