

Product Datasheet

Mouse Myonectin ELISA Kit (orb1212388)

Description	Mouse Myonectin, CTRP15 ELISA Kit
Reactivity	Mouse
Range	0.04-15ng/ml
Tested Applications	ELISA
Target	ERFE
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	36,265
Uniprot ID	Q6PGN1
Solubility (25°C)	16ng/ml
Sensitivity	0.02ng/ml
Entrez	227358
Expiration Date	Please enquire.