

Product Datasheet

Mouse Neuromedin U ELISA Kit (orb1212392)

Description	Mouse Neuromedin U, NMUELISA Kit
Reactivity	Mouse
Range	0.5-100ng/mL
Tested Applications	ELISA
Target	NMU
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	19429
Uniprot ID	Q9QXK8
Solubility (25°C)	128ng/mL
Sensitivity	0.224ng/mL
Entrez	56183
Expiration Date	Please enquire.