

Product Datasheet

Mouse Fibrillin-1 ELISA Kit (orb1212397)

Description	Mouse Fibrillin-1, FBN1 ELISA Kit
Reactivity	Mouse
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	FBLN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	78,033
Uniprot ID	Q08879
Solubility (25°C)	64ng/ml
Sensitivity	0.013ng/ml
Entrez	14114
Expiration Date	Please enquire.