

Product Datasheet

Mouse Annexin A1 ELISA Kit (orb1212400)

Description	Mouse Annexin A1, ANXA1 ELISA Kit
Reactivity	Mouse
Range	0.2-70ng/mL
Tested Applications	ELISA
Target	ANXA1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	38,734
Uniprot ID	P10107
Solubility (25°C)	80ng/mL
Sensitivity	0.118ng/mL
Entrez	16952
Expiration Date	Please enquire.