

Product Datasheet

Mouse S100 Calcium Binding Protein A8/Calgranulin A ELISA Kit (orb1212416)

Description	Mouse S100 Calcium Binding Protein A8/Calgranulin A, S100A8 ELISA Kit
Reactivity	Mouse
Range	0.2-70ng/mL
Tested Applications	ELISA
Target	S100A8
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	10,295
Uniprot ID	P27005
Solubility (25°C)	80ng/mL
Sensitivity	0.121ng/mL
Entrez	20201
Expiration Date	Please enquire.