

Product Datasheet

Mouse Fructosamine-3-kinase ELISA Kit (orb1212417)

Description	Mouse Fructosamine-3-kinase, FN3K ELISA Kit
Reactivity	Mouse
Range	30-9000ng/L
Tested Applications	ELISA
Target	FN3K
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	35,032
Uniprot ID	Q9ER35
Solubility (25°C)	9600ng/L
Sensitivity	14.36ng/L
Entrez	63828
Expiration Date	Please enquire.