

Product Datasheet

Mouse Acetylcholinesterase ELISA Kit (orb1212418)

Description	Mouse Acetylcholinesterase, ACHE ELISA Kit
Reactivity	Mouse
Range	0.5-150ng/ml
Tested Applications	ELISA
Target	ACHE
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	68,169
Uniprot ID	P21836
Solubility (25°C)	160ng/ml
Sensitivity	0.27ng/ml
Entrez	11423
Expiration Date	Please enquire.