

Product Datasheet

Mouse Lectin-like Oxidized Low-density Lipoprotein Receptor 1 ELISA Kit (orb1212419)

Description	Mouse Lectin-like Oxidized Low-density Lipoprotein Receptor 1, OLR1 ELISA Kit
Reactivity	Mouse
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	OLR1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	41,643
Uniprot ID	Q9EQ09
Solubility (25°C)	32ng/ml
Sensitivity	0.023ng/ml
Entrez	108078
Expiration Date	Please enquire.