

Product Datasheet

Mouse Kallikrein 8 ELISA Kit (orb1212422)

Description	Mouse Kallikrein 8, KLK8 ELISA Kit
Reactivity	Mouse
Range	20-6000ng/L
Tested Applications	ELISA
Target	KLK8
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	28,524
Uniprot ID	Q61955
Solubility (25°C)	6400ng/L
Sensitivity	10.33ng/L
Entrez	259277
Expiration Date	Please enquire.