

Product Datasheet

Mouse Pulmonary Surfactant-associated Protein D ELISA Kit (orb1212424)

Description	Mouse Pulmonary Surfactant-associated Protein D, SP-D ELISA Kit
Reactivity	Mouse
Range	0.5-200ng/ml
Tested Applications	ELISA
Target	SFTPD
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	37,688
Uniprot ID	P50404
Solubility (25°C)	240ng/ml
Sensitivity	0.3ng/ml
Entrez	20390
Expiration Date	Please enquire.