

Product Datasheet

Mouse Glycogen Synthase Kinase-3 Beta ELISA Kit (orb1212425)

Description	Mouse Glycogen Synthase Kinase-3 Beta, GSK3B ELISA Kit
Reactivity	Mouse
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	GSK3B
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	46,710
Uniprot ID	Q9WV60
Solubility (25°C)	40ng/ml
Sensitivity	0.024ng/ml
Entrez	56637
Expiration Date	Please enquire.