

Product Datasheet

Mouse Endothelial Cell-specific Molecule 1 ELISA Kit (orb1212430)

Description	Mouse Endothelial Cell-specific Molecule 1, ESM-1 ELISA Kit
Reactivity	Mouse
Range	5-2000ng/L
Tested Applications	ELISA
Target	ESM1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,043
Uniprot ID	Q9QYY7
Solubility (25°C)	2400ng/L
Sensitivity	2.53ng/L
Entrez	71690
Expiration Date	Please enquire.