

Product Datasheet

Mouse Beta-defensin 4 ELISA Kit (orb1212436)

Description	Mouse Beta-defensin 4, DEFB4 ELISA Kit
Reactivity	Mouse
Range	0.05-10ng/ml
Tested Applications	ELISA
Target	DEFB4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	7,129
Uniprot ID	P82019
Solubility (25°C)	12.8ng/ml
Sensitivity	0.021ng/ml
Entrez	56519
Expiration Date	Please enquire.