

Product Datasheet

Mouse Ficolin 2 ELISA Kit (orb1212438)

Description	Mouse Ficolin 2, FCN2 ELISA Kit
Reactivity	Mouse
Range	70-15000ng/L
Tested Applications	ELISA
Target	FCN2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	33,984
Uniprot ID	070497
Solubility (25°C)	16000ng/L
Sensitivity	34.65ng/L
Entrez	14134
Expiration Date	Please enquire.