

Product Datasheet

Mouse Glucagon-like Peptide 1 Receptor ELISA Kit (orb1212448)

Description	Mouse Glucagon-like Peptide 1 Receptor, GLP1R ELISA Kit
Reactivity	Mouse
Range	15-3000ng/L
Tested Applications	ELISA
Target	GLP1R
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	53,028
Uniprot ID	035659
Solubility (25°C)	3200ng/L
Sensitivity	7.05ng/L
Entrez	14652
Expiration Date	Please enquire.